

Catalog



Dental press furnace

NanYang LianDong Biotechnology Co.,Ltd.

Email: admin@nyliandong.com
Web: www.xianghuism.com

Tel: +86-0377-66077008

Brief Introduction

LD-P-V1200 Press Furnace for press unit is a kind of electric heating press, which can makes the press force more even. Besides, it can accept the soft investment material, like Ivoclar SPEED, so the dental veneer can get much thinner, more better transparency, better sealness.



NanYang LianDong Biotechnology Co.,Ltd.

Email: admin@nyliandong.com
Web: www.xianghuism.com

Tel: +86-0377-66077008

Product Features

- Easy to operate, compact size, one of the world's smallest sizes porcelain furnace at present.**
- Automatically high-quality temperature adjustment function ensures that the actual temperature error is within $\pm 1^{\circ}$.**
- Precision stepper motor drive, smooth free operation and no jitter functioning, noise-free, speed programmable**
- Spiral quartz chamber of furnace.**
- The tray has no noise and its speed is programmable.**
- USB port facilitates connection with computer service programs.**
- True color touch screen, easy to operate.**
- 99 built-in programs, and customized programs**
- Optical status display indication.**

NanYang LianDong Biotechnology Co.,Ltd.

Email: admin@nyliandong.com
Web: www.xianghuism.com

Tel: +86-0377-66077008

Specification

Model	LD-P-V1200
Voltage	110V/220C,50/60HZ
Power	Max.1.5KW
Electric current	14.5ampere@110V
Max. Temperature	1200°C
Working Temperature	≤1150°C
Heating rate	100°C/min max.
Vacuum degree	-98Kpa
Thermocouple	K type
Heating element	Fe-Cr-Al-Mo wire
Packing size	41*54*85cm(Press furnace);41*32*28cm(Vacuum Pump)
Packing Weight	27kg(Press Furnace),10kg(Vacuum Pump)

NanYang LianDong Biotechnology Co.,Ltd.

Email:admin@nyliandong.com
Web:www.xianghuism.com

Tel:+86-0377-66077008


Accessories

Items	Quantity
Vacuum Pump	1 set
Embedding Ring	2 pcs(100g/200g)
Press Rod	2pcs
Fire Pan	1pc
Cooling Rack	1pc
Tong	1pc

NanYang LianDong Biotechnology Co.,Ltd.

Email:admin@nyliandong.com
Web:www.xianghuism.com

Tel:+86-0377-66077008



Application Case

NanYang LianDong Biotechnology Co.,Ltd.

Email: admin@nyliandong.com
Web: www.xianghuism.com

Tel: +86-0377-66077008



NanYang LianDong Biotechnology Co.,Ltd.

Email: admin@nyliandong.com
Web: www.xianghuism.com

Tel: +86-0377-66077008



NanYang LianDong Biotechnology Co.,Ltd.

Email: admin@nyliandong.com
Web: www.xianghuism.com

Tel: +86-0377-66077008

Products show



NanYang LianDong Biotechnology Co.,Ltd.

Email: admin@nyliandong.com
Web: www.xianghuism.com

Tel: +86-0377-66077008